



Ion/Molecule Reactions in MS

Permanent dipole and polarizability of atoms and molecules




Atom, molecule	Permanent dipole, μ (Debye)	Polarizability, α ($\text{\AA}^3$)
Rare gas atoms		
He	0	0.21
Ne	0	0.41
Ar	0	1.64
Kr	0	2.48
Xe	0	4.02
Diatomic		
H ₂	0	0.81
D ₂	0	0.8
O ₂	0	1.57
N ₂	0	1.74
CO	0.1	1.94
NO	0.16	1.7
HCl	1.03	2.58
HBr	0.79	3.49
Triatomic molecule		
CO ₂	0	2.59
CS ₂	0	8.08
H ₂ O	1.84	1.45
N ₂ O	0.17	2.92
H ₂ S	1.02	3.61
SO ₂	1.62	3.78
Polyatomic molecule		
CH ₄	0	2.56
CCl ₄	0	10.24
C ₂ H ₄	0	4.1
C ₂ H ₆	0	4.39
n-C ₄ H ₁₀	0	8
i-C ₄ H ₁₀	0	8
C ₆ H ₆	0	9.99
NH ₃	1.47	2.16
H ₂ CO	2.31	2.81
CH ₃ OH	1.7	3.25
(CH ₃) ₂ CO	2.85	6.11



The measure of “softness” of an electron cloud is quantitatively given by the “polarizability” (α).

离子靠近极性分子的过程



Types of Ion/molecule reactions


- Electron Transfer Reactions
- Proton Transfer Reactions
- Hydride Transfer Reactions
- Clustering Reactions
- Negative Ion Formation
- SN2 displacement reactions

Electron Transfer Reactions



Proton Transfer Reactions

MALDI (Matrix-Assisted Laser Desorption/Ionization, 基质辅助激光解吸电离)



<https://www.creative-proteomics.com/>

1. 样品制备

选择基质：选择适合分析物的基质，常用基质包括 α -氰基-4-羟基肉桂酸 (CHCA) 和二羟基苯甲酸 (DHB)。

混合样品与基质：将样品与基质溶液混合，通常比例为1:1000到1:10000。

点样：将混合液滴在MALDI靶板上，待溶剂挥发后形成共结晶。

2. 离子源准备

靶板装载：将点样后的靶板放入MALDI离子源。

真空系统：启动真空系统，确保离子源内部达到所需真空度。

3. 激光照射

激光参数设置：根据样品特性调整激光波长（通常为337 nm）和能量。

激光照射：激光脉冲照射样品-基质共结晶区域，使基质吸收能量并解吸样品分子。


4. 离子化过程

能量转移：基质吸收激光能量并转移给样品分子，使其离子化。


5. 离子生成

样品分子形成正离子或负离子，通常为单电荷离子。

化学电离 (CI) 的质子转移:




Hydride Transfer Reactions



Negative Ion Formation

Electron capture:





Ion pair formation:



Charge transfer:



Electron capture reactions frequently used for CI



Electron affinity (eV)

Chemical species	Electron affinity
H	0.75
O	1.46
F	3.4
Cl	3.62
Br	3.37
I	3.07
NO	0.024
NO ₂	2.43
NO ₃	3.7
N ₂ O	0.22
O ₂	0.5
O ₃	2.05
OH	1.83
SCN	2.17
SF ₆	0.6
SF ₅	2.8
7,7',8,8'-tetracyanoquinodimethane (TCNQ)	2.83
Nitrobenzene	0.98
Tetracene	0.88
Azulene	0.66
Anthracene	0.55
Naphthalene	0.14

A typical example of SN2 reaction is the halide ion displacement reaction:

